


## CGQ BioR


Online Biomass Monitoring for Bioreactors



The CGQ BioR is a sensor for online biomass monitoring in various types and scales of bioreactors.


### The CGQ BioR

---



Backscatter measurements are used to monitor the biomass non-invasively through the wall of the bioreactor.


### The CGQ BioR Sensor & Its Principle of Measurement





The CGQ BioR is mounted to the outside of the reactor vessel and is compatible with the majority of bioreactors (single wall and double wall).

### Exemplary Pictures of CGQ BioR Sensors Mounted to Different Types of Bioreactors



The CGQ BioR creates significant value for its users by saving time and costs as well as creating detailed understanding of the bioprocess.

## CGQ BioR Features & Benefits



### Features & Benefits

#### **Non-Invasive Plug & Play Technology**

*No cleaning or autoclaving, no ports blocked, quickly install/uninstall at any time*

#### **Compatible with Various Bioreactors**

*Fully compatible with most bioreactors from all common manufacturers*

#### **Detailed Microbial Growth Kinetics in Real-Time**

*Real-time analysis of important process parameters in highest precision*

#### **One Sensor for All Needs**

*Covers the OD range from 0.5 to 300 in two measurement modes*

#### **Saves Time & Costs**

*No manual sampling for biomass quantification needed, increases productivity and reduces personnel costs*

#### **Monitoring at All Times**

*Continuous measurements: during the day, at night and over the weekend*



The CGQ BioR has several clear advantages over currently used invasive biomass sensor for bioreactors.


### Comparison: CGQ BioR vs. Invasive Biomass Sensor

	CGQ BioR	Invasive Biomass Sensor
Cleaning & Autoclaving	<ul style="list-style-type: none"><li>– Non-invasive sensor that <b>does not need to be autoclaved or cleaned</b></li></ul>	<ul style="list-style-type: none"><li>– <b>Needs to be cleaned and autoclaved</b> with the vessel after every use</li></ul>
Available Ports	<ul style="list-style-type: none"><li>– <b>No ports blocked</b> since the BioR is attached to the outside of the glass vessel</li></ul>	<ul style="list-style-type: none"><li>– <b>Requires a port</b></li></ul>
Flexibility	<ul style="list-style-type: none"><li>– <b>Can be installed/uninstalled at any given time</b> during the fermentation (<b>plug and play</b>)</li></ul>	<ul style="list-style-type: none"><li>– <b>Must be installed before the experiment</b> is started</li></ul>
Vessel Compatibility	<ul style="list-style-type: none"><li>– <b>Compatible with most vessel types and sizes</b></li><li>– Simply attach to the glass wall or a glass window</li></ul>	<ul style="list-style-type: none"><li>– Can often only be used <b>for one vessel size</b> (limited by probe length)</li></ul>
OD Range	<ul style="list-style-type: none"><li>– Standard Mode (521 nm): OD 0.5-50*</li><li>– High Cell Density Mode (940 nm): OD 15-300*</li></ul>	<ul style="list-style-type: none"><li>– Depending on manufacturer</li><li>– Usually <b>limited to a specific OD range</b></li></ul>

\*Depending on vessel type/size, media, organism and other factors

The CGQ BioR shows good correlation with offline biomass data such as OD<sub>600</sub> and Cell Dry Weight.

CGQ BioR & Offline Biomass Measurements (1/2): *Corynebacterium glutamicum*




*Corynebacterium glutamicum*, CGXII Medium, DASGIP Bioblock (1,8 l), 30 °C



The CGQ BioR shows good correlation with offline biomass data such as Cell Count and Biovolume.

CGQ BioR & Offline Biomass Measurements (2/2): *Corynebacterium glutamicum*




*Corynebacterium glutamicum*, CGXII Medium, DASGIP Bioblock (1,8 l), 30 °C





Using a calibration file, the CGQ BioR is able to directly convert backscatter measurements to OD values.


### CGQ BioR and Offline OD Measurements: *Bacillus subtilis*



*Bacillus subtilis*, Minimal Medium, New Brunswick Glass Bioreactor (3 l), 37 °C

The CGQ BioR is ideal for screening experiments in bioreactors such as strain or media comparisons.


### CGQ BioR Measurement: Different *Escherichia coli* Strains



*Escherichia coli*, HCDC Medium, Infors Labfors (7,5 l), 37 °C

The CGQ BioR allows for non-invasive online biomass monitoring in glass bioreactors with a double wall.


### CGQ BioR Measurement: *Escherichia coli*



*Escherichia coli* BL21, M9 Medium, Applikon Glass Bioreactor (7,5 l), Room Temperature

The high data density of CGQ BioR measurements enables the user to detect and visualize process events like metabolic shifts.

### CGQ BioR Measurement: *Saccharomyces cerevisiae* (Diauxic Growth)




*Saccharomyces cerevisiae*, YPD Medium, Applikon Glass Bioreactor (7,5 l), Room Temperature



The CGQ BioR allows the user to closely follow the growth but also the dying of his/her cultures.

CGQ BioR Measurement: *Amycolatopsis* sp. ATCC 39116



*Amycolatopsis* sp. ATCC 39116, Glucose Yeast & Malt Medium, Sartorius BIOSTAT b plus (2 l), 45 °C